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Thank you for purchasing a HORIBA ISFET pH electrode. This product is less fragile and more easy to handle than existing glass electrodes because of semiconductor, and can be used with pH meters for glass electrodes (D-10, D-20, D-50, F-10, F-20, F-50 series). Read this manual carefully before use.

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■ Warranty and Responsibility

The product (the converter excluding batteries, sensor, and accessories) shall be repaired or replaced free of charge should any malfunction occur during the course of normal use if returned within the warranty period (12 months from date of purchase). Return the converter with proof of purchase date.

This warranty does not cover damages caused by accidents, disasters, misuse, carelessness, and unauthorized disassemble and remodeling.

The sole and exclusive remedy of HORIBA under this limited warranty shall be repair or replacement of the product within a commercially reasonable time from notification.

In any case, HORIBA shall not be responsible for customer's passive damages, or damages which may arise from claims instituted by a third party against customers.

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1. Cautions on use

- Handle carefully the sharp electrode tip to not get hurt.
- If calibration standard solution gets on the skin, wash off it completely. In case that the solution gets into the eyes, continue to flush the eyes with running water for at least 15 minutes. And after washing the eyes thoroughly, get medical help immediately.
- Do not eat the whole, or the part around where the electrode sensor touched, of measured samples. Some of the sensor materials may generate environmental endocrine disrupters.
- Do not disassemble the product. This may cause malfunction.
- Do not drop the electrode. And do not apply excessive force to the sensor and converter.
- Do not scratch or apply excessive force to the sensing part and liquid junction.
- Do not touch the sensing part and sensor terminal. This may cause malfunction.
- Do not allow the electrode connector and sensor terminal to come in contact with water or wet hands.
- Before replacing the sensor or batteries, turn off the power and wipe water off the electrode.
- For the first use, or use after a long interval (one week or longer), immerse the sensor tip in 3.33 mol/L KCl solution (#300) for 30 minutes or longer before use.
- Make sure to use the electrode at a temperature within the specified usable range (sensor: 0 °C to 60 °C , converter: 0 °C to 40 °C).
- The product should not be subject to rapid temperature change.
- Keep the sensor away from direct sunlight or strong light during calibration and measurement, as accuracy will be affected adversely.
- Avoid continuous use in acids or alkali sample of 0.01 mol/L or thicker. Failure to comply may cause decline in the performance and life period of the electrode.
- When sticking the sensor into a solid sample, take care not to apply any bending or twisting force to the sensor. This may cause malfunction.
- Dirt on the sensing part or liquid junction will affect accuracy adversely. Keep the sensor clean referring to "7. Maintenance/ Cleansing the sensing part" (page 11).
- Store the product in dry conditions. Storage area should be kept away from high temperature, high humidity, and direct sunlight.
- Do not use organic solvents such as acetone, tetrahydrofuran (THF) for cleaning. The sensor materials include resins, including ABS resin, epoxy resin, and polyethylene resin.
- White powder or solution may be generated at the liquid junction. This is leaked internal solution of the reference electrode, and not a symptom of trouble. Rinse the sensor with ion-exchange water before use.
- The sensor is a consumable. If it is damaged or deteriorated, replace a new sensor. (The sensor cannot be repaired.)
- Make sure to detach the sensor before battery replacement. Failure to comply may cause malfunction.
- Misusing batteries may cause an internal solution leak or a damage of the batteries. Follow the instruction below.
 - CR2032 batteries are not rechargeable.
 - Remove the batteries from the converter before prolonged nonuse.
 - If the battery solution leaks, wipe away the spilled solution from the battery case before replacing the batteries.
- For child's safety, keep batteries away from children. If your child does swallow a battery, take him or her to the emergency room immediately.
- If the Low Batt. LED blinks, the battery voltage has been low. In this case, replace the batteries as soon as possible. When the Meas. LED is turned off, electrode performance is not guaranteed.
- Do not throw batteries into fire.


2. Specifications and Components

■ Packaged Contents


- Sensor (0131): 1
- 3.33 mol/L KCl solution #300: 1
- Converter: 1
- Instruction manual (this book): 1
- Battery CR2032: 2
- Screw driver: 1

■ Components

● Sensor (0131)

● Converter

■ Specifications


Electrode model: 0030-10D
 Sensor model: 0131
 Measurement range: pH 0 to 14
 Usable temperature range: [Sensor] 0 °C to 60 °C , [Converter] 0 °C to 40 °C
 Storage temperature range: 0 °C to 40 °C
 Electrode length: 190 mm (excluding cable)
 Sensing device: ISFET (Ion-sensitive field-effect transistor)
 Internal electrode: Silver/ silver chloride electrode
 Reference electrode internal solution: Carboxymethyl cellulose (CMC), glycerin, supersaturated KCl
 Liquid junction material: porous sintered polyethylene
 Sensor material: Tantalum pentoxide, platinum, acrylonitrile-butadiene-styrene (ABS) resin, epoxy resin, polyethylene resin
 Height of sensing part: 16 mm (from the tip of electrode)
 Application: D-10, D-20, D-50, F-10, F-20, F-50 series
 Power: 6 V DC CR2032 × 2, Auto power OFF
 LED indicator: Meas. LED (green blink while the power is ON), Low Batt. LED (red blink when the battery voltage is low)
 Cable length: 1 m

3. Preparation


For the first use, prepare the ISFET pH electrode as follows:

■ Setting batteries


- 1** Unhook and draw the converter cover slowly.

- 2** Set two new batteries (CR2032) by slide them under the positive contacts with the positive poles of the batteries up.

- 3** Attach and hook the converter cover to the converter by making the projection (A) fit in the hole (a) on the hook.


■ Assembling sensor and converter


Note

Do not touch the sensor terminal. It may be damaged by static electricity.


- 1** Remove the end cap from the sensor, and the seal from the liquid junction. Make sure that the O-ring is attached to the sensor.


- 2** Turn the nut to the left to remove it from the converter.

- 3** Insert the sensor to the converter slowly checking that the O-ring is untwisted, and match the Δ marks on the sensor and converter. Then turn the nut to the right to fix the sensor.

- 4** Immerse the sensing part and liquid junction completely into 3.33 mol/L KCl solution (#300) for 30 minutes or longer.

Before measurement or calibration, prepare the ISFET pH electrode as follows:

■ Connecting to pH meter



- 1** Insert the electrode connector into the connector port on the meter, after aligning with the pin. Do not insert the connector unless it is aligned properly with the port.
- 2** Press the electrode connector into the port on the meter, while turning the connector to the right.
- 3** Slide the connector cover over the electrode connector. Then, push the cover in straight until it comes in light contact with the meter case. Do not turn the cover.
- 4** Insert the temperature connector securely into the jack on the meter.

4. Calibration

Note

ISFET needs the zero adjustment procedures before calibration with standard solution.



- 1 Immerse the sensing part and liquid junction completely into 3.33 mol/L KCl solution (#300) for 10 minutes or longer. After that, rise them with pure water and wipe them with filter paper.
- 2 Press the power switch to turn ON the power. The Meas. LED blinks green.

Note

Automatic power OFF works in approximately 1 hour after power ON.

- 3 Keep the electrode in pH7 standard solution, with the sensing part and liquid junction completely immersed. If necessary, swing the sensor gently in the solution to remove bubbles near the sensing part and liquid junction.
- 4 Adjust the zero point of the ISFET pH electrode by turn the zero adjustment volume with a slot driver, to keep pH meter readings within the range shown right.

In mV measurement mode:
 0 ± 10 mV
In pH measurement mode:
 ± 0.2 pH for the pH value of pH7 standard solution

Note

Before adjustment in pH measurement mode, clear the existing calibration values and check that "CLR" is displayed on the screen. See the instruction manual of the pH meter for the details.



- 5 Calibrate the ISFET pH electrode according to the instruction manual of the pH meter you use.

5. Measurement



- 1 Wash the sensing part and liquid junction with pure water and dry them with filter paper, prior to every immersion into standard/sample solution.
- 2 Press the power switch to turn ON the power. The Meas. LED blinks green.

Note

Automatic power OFF works in approximately 1 hour after power ON.

- 3 For liquid samples:
Immerse the sensing part and liquid junction to the sample completely. Stir the electrode lightly to remove bubbles near these parts, if necessary.
For solid samples:
Press the electrode tip into the sample slowly. Make sure to keep the sensing part close to the sample without a gap.
- 4 Measure the sample according to the instruction manual of the pH meter you use.

Note

Do not eat the whole, or the part around where the electrode sensor touched, of measured samples. Some of the sensor materials may generate environmental endocrine disrupters.

● Recommended conditions for accurate measurement

For accurate measurement, note the items below.


- Stir the sample gently with a stirrer during measurement or calibration.
- Standard solutions used for calibration and measurement samples should be the same in temperature.
- Rinse the electrode with the sample before measurement.
- Calibrate the ISFET pH electrode before every measurement.
- Before use, immerse the sensor into 3.33 mol/L KCl solution (#300) for 30 minutes or longer.

6. Storage

- 1 Press the power switch to turn OFF the power. The Meas. LED turns off.
- 2 Wash the electrode well with pure water to remove sample completely and wipe it with filter paper, and store it in dry conditions.

Note

Avoid storing the electrode in a hot place or a highly humid place.
Keep the electrode dry and out of direct sunlight.

7. Maintenance


■ Cleansing the sensing part

Dirt on the sensing part may affect the electrode response, sensitivity, or a measurement precision. If the electrode is very dirty and cannot be washed clean using pure water, clean it as follows.

- 1 Clean the sensing part by wiping it with cotton gauze containing cleansing agent.
- 2 Rinse the electrode with pure water completely.

Note

Be sure to calibrate the electrode after cleansing the sensing part.

Type of dirt	Cleansing agent
General	Diluted neutral cleansing agent
Oily	Diluted neutral cleansing agent, or alcohol
Inorganic (AgCl, etc.)	Cleansing agent #220

Note


Do not use organic solvent (eg. acetone, THF), strong acid or alkali. Using these chemicals may damage the electrode.

■ Replacing the sensor

If any of the following symptoms is found after reponsive-part cleaning, replace the sensor in the procedures below.

- Sensitivity or response speeds do not recover.
- The electrode can not be adjusted in the way mentioned in "4. Calibration" (page 10).
- The electrode can not be calibrated with standard solutions.

- 1 Turn OFF the power.
- 2 Turn the nut to the left to remove it. And pull out the sensor from the converter.

- 3 Attach a new sensor referring to "3. Preparation/ Assembling sensor and converter" (page 9).


Note

Do not touch the sensor terminal.

■ Replacing batteries

When battery voltage is low, the Low Batt. LED blinks red. (The LED is off during battery shutoff.) In this case, replace batteries as follows:

- 1 Turn OFF the power.
- 2 Turn the nut to the left to remove it. And pull out the sensor from the converter.

- 3 Replace batteries with new ones referring to "3. Preparation/ Setting batteries" (page 9). If used batteries are set, remove it by using a ball-point pen or the like.

- 4 Attach the sensor referring to "3. Preparation/ Assembling sensor and converter" (page 9).

Note

- Replace two batteries with new ones (CR2032) at the same time.
- Do not touch the sensor terminal.
- Make sure to remove the sensor prior to the converter cover. Removing the cover with the sensor attached causes malfunction.

For more information, please contact us:

[ExpotechUSA](#)
[10700 Rockley Road](#)
[Houston, Texas 77099](#)
[USA](#)

[281-496-0900 \[voice\]](#)

[281-496-0400 \[fax\]](#)

E-mail: sales@expotechusa.com

Website: www.ExpotechUSA.com