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## Operation Manual for X-52 Simulator

Carefully read this manual before using the X-52 Simulator. X-52 Simulator gives theoretical equivalent input of electrodes for Conductivity measuring instruments and checks linearity and temperature compensation functions of the instrument. The simulator makes equivalent resistance output for conductivity and thermistor equivalent temperature resistance output for temperature.

### HORIBA's Warranty and Responsibility

X-52 Simulator is covered by HORIBA's warranty for a period of one (1) year, under normal use. Although unlikely, if any trouble attributable to HORIBA should occur during this period, necessary exchange or repairs shall be conducted by HORIBA, free of charge. The warranty does not cover the following:

- Any trouble or damage attributable to actions or conditions specifically mentioned to be avoided in the operation manuals
- Any trouble or damage attributable to use of the pH meter in ways or for purposes other than those described in the operation manuals
- Replacement of consumables and accessories, or use of parts not specified by HORIBA
- Any alteration to the external appearance of this Simulator attributable to scratches, dirt, etc. occurring through normal use
- Any trouble or damage attributable to repair or modification by any party not authorized by HORIBA

HORIBA also shall not be liable for any damages resulting from any malfunctions of this product, any erasure of data, or any other uses of this product.

### CE Marking

This equipment is in conformity with the following directives and standards:

Directives: The EMC Directives 89/336/EEC  
The Electrical Product Safety Directive 73/23/EEC

Standards: EN61326: 1997+A1: 1998  
(EMISSION: Class B, IMMUNITY Category: Minimum Requirement)  
EN61010-1: 1993

### Installation Environment

This product is designed for the following environment.

- Pollution degree 2
- Measurement categorie I

WARNING: Do Not use the equipment for measurements within measurement categories II, III and IV.

### FCC Warning

This equipment has been tested and found to comply with the limits for a Class A digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment.

This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference in which case the user will be required to correct the interference at his own expense.

### Safety precautions for use

- Do not give physical shock to the Simulator like dropping or hitting.
- Do not immerse the Simulator into alcohol, organic solvent, strong acid, strong alkaline, and other similar solutions. The simulator contains ABS resin, acrylic resin and various rubber products in its body.
- Do not use a hair-dryer for drying the Simulator. When the Simulator is wet, wipe it using soft cloth.
- Perform the key operation by the fingers, not by the hard object like metal stick or rod.
- To disconnect the cable, pull them out with holding the connector part. Do not pull the cable part; it may cause a breakage.
- Do not remove the battery gasket or twist it.
- When opening the battery case, make sure that no foreign matter is attached to the battery gasket.
- Do not use any unspecified batteries; it may cause a breakage.
- Be sure to use packaged cable. The accuracy is guaranteed by the combination of Simulator and cable.

### Auto power OFF function

This Simulator automatically turns off, if it is left without any key operation. To turn the power OFF manually, press ON/OFF key for one second or more.

Termination of output	After one minute without key operation
Power OFF	After five minutes without key operation

### Storage and operating environment

- Storage temperature: 0°C to 45°C,  
Operating temperature: 5°C to 40°C
- Place the Simulator where the relative humidity is kept below 80% and no dew is condensed.
- Avoid operating the Simulator in the following kind of locations.
  - Dusty locations
  - Location with severe vibration
  - Location with direct sunlight

- Location exposed to corrosive gas
- Location near to an air-conditioner or heater
- Location exposed to direct wind

### Transportation of the X-52

- To transport the Simulator, use the packaging box at the delivery. Transportation by any unspecified packing methods may cause a breakage.


### Disposition

- When you dispose the Simulator, be sure to dispose it as an industrial waste.

## 1 Package contents


- Simulator (X-52) 1 pc
- Cable (1 m) 1 pc
- Dry-cell battery (AA alkaline) 2 pcs
- Strap 1 pc
- Soft case 1 pc
- Operation manual 1 sheet

## 2 Part names


Keys	Description
	ON/OFF key Turns ON/OFF the power
	OUTPUT key Outputs simulation input
	MODE key Toggles output mode
	UP key Toggles simulation input value
	DOWN key Toggles simulation input value

## 3 Procedures for replacing battery

## 4 Connection

Connect the end of connector of cable to the desired output portion, and other end to the corresponding electrode connection port of instrument.


## 5 Check

Turn on the power of Simulator with ON/OFF key.

- When certifying the instrument, keep the room temperature in the range of  $\pm 5^\circ\text{C}$  of the temperature in X-52 delivery inspection.  
\*: Temperature in X-52 delivery inspection is described in the inspection report.
- It is necessary to keep the Simulator placed in environment of  $\pm 5^\circ\text{C}$  of the temperature in X-52 delivery inspection for one hour to secure the output accuracy.
- Follow the "Operational Qualification Confirmation Report", when performing OQ (Operational Qualification).

**Note**▶▶ When the battery is running down, "Err" is displayed, and the Simulator does not start.


### Caution for the operation during outputting simulation value

- If OUTPUT key is pressed or mode is changed with MODE key, outputting simulation input value stops.
- Even if simulation input value is toggled, outputting simulation input value continues unless output mode is changed.


### 5.1 Conductivity check

1. Set the cell constant to  $1.000 \times 100 \text{ m}^{-1}$  on the instrument.
2. Select COND Span with MODE key on the Simulator.
3. Start outputting simulation input value with OUTPUT key.
4. Input span simulation value in the order of 19.00 S/m, 190.0 S/m, 19.00 mS/m, and 1.900 mS/m to the Simulator, then check the value on the instrument.

- 2. Select COND Span.
- 4.-① Output simulation input value.
- 4.-② Toggle simulation input value.


5. Select COND Linear with MODE key on the Simulator.
6. Start outputting simulation input value with OUTPUT key.
7. Input linear simulation value in the order of 10.00 S/m, 1.000 S/m, 100.0 S/m, 10.00 mS/m, 1.000 mS/m, 100.0  $\mu\text{S/m}$ , 10.00  $\mu\text{S/m}$ , and 0.000 mS/m to the Simulator, then check the value on the instrument.

- 5. Select COND Linear.
- 7.-① Output simulation input value.
- 7.-② Toggle simulation input value.

### 5.2 Temperature display check

1. Select Temp Set with MODE key on the Simulator.
2. Start outputting simulation input value with OUTPUT key.
3. Input temperature simulation value in the order of 0.0°C, 25.0°C, 30.0°C, 60.0°C and 100.0°C to the Simulator, then check the value on the instrument.

- 1. Select Temp Set.
- 3.-① Output simulation input value.
- 3.-② Toggle simulation input value.

## 6 Specifications

Output parameter	Output value	Accuracy (Varies depending on temperature used.)		
		$\pm 5^\circ\text{C}$ of the temperature in X-52 delivery inspection(*)	$5^\circ\text{C}$ to $40^\circ\text{C}$	
Output	COND Span	19.00 S/m	$\pm 0.1 \text{ S/m}$	$\pm 0.8 \text{ S/m}$ $\pm 1 \text{ digit}$
		1.900 S/m to 1.900 mS/m	$\pm 0.1\% \pm 1 \text{ digit}$	$\pm 0.2\% \pm 1 \text{ digit}$
	COND Linear	10.00 S/m	$\pm 0.1 \text{ S/m}$	$\pm 0.4 \text{ S/m}$ $\pm 1 \text{ digit}$
		1.000 S/m to 10.00 $\mu\text{S/m}$	$\pm 0.1\% \pm 1 \text{ digit}$	$\pm 0.2\% \pm 1 \text{ digit}$
		0.000 mS/m	-	-
	Temperature	0.0/ 25.0/ 30.0/ 60.0/ 100.0°C	$\pm 0.1^\circ\text{C}$	$\pm 0.2^\circ\text{C}$
Ambient temp.	5°C to 40°C			
Ambient humidity	0% to 80% (No dew condensation)			
Power	Two AA alkaline dry-cell batteries (Approx. 20 hours with continuous operation)			
Dimensions	170(H) × 89(W) × 40(D) mm			

(\*): Temperature in X-52 delivery inspection is described in the inspection report.

**Note**▶▶ This Simulator supplies simulation input specific to HORIBA Conductivity meter. Therefore, the simulator output cannot be measured by general-purpose digital multi-meter.

For more information, please contact us:

ExpotechUSA

10700 Rockley Road  
Houston, Texas 77099  
USA

281-496-0900 [voice]

281-496-0400 [fax]

E-mail: sales@expotechusa.com

Website: www.ExpotechUSA.com