


HORIBA

OIL CONTENT ANALYZER OCMA-350





Easy-to-read backlit LCD

The backlit LCD is easy to read under any conditions, including bright sunlight. The digital display eliminates individual reading differences.

Three measurement units

The OCMA-350 measures in mg/l for oil in water, mg/kg for oil in soil, and Abs. for the absorbency of oil, so you can always choose the proper unit for your sample measurement.

Cell control

To introduce your sample into the removable cell, just pull the cell out of its holder. It's easy to avoid leaks and spills when you have this much control. The removable cell is also easy to maintain.



**Plug it in, and it's on stand-by.
The press of just one button
starts the measuring process.
You can even use your PC to
simplify operation.**


Pushbutton simplicity, perfect for field analysis.

Compact and simple, HORIBA's OCMA-350 oil content analyzer delivers quick, accurate measurements in the lab or in the field. It efficiently measures the level of oil in water or soil for environmental applications, checks for residual oil on semiconductor parts which have been cleaned, and measures oil on any industrial surface.

Using the OCMA-350 is almost effortless. Simply inject the extract from your sample, dissolved in Horiba's S-316 solvent, into the measuring cell, and press one button. You can even use your PC to simplify operation.

Recommended Applications

- Environment: Surveying water quality and hazardous waste sites
- Industry: Monitoring wastewater influent and effluent
- Marine transportation: Checking bilge and ballast discharge
- Petroleum processing: Monitoring the efficiency of oil/water separation processes
- Oil depots: Monitoring the discharge produced when cleaning storage tanks
- Automotive: Monitoring waste-water discharge from service stations
- Quality Assurance: Measuring residual oil on textiles and metal parts.



All the right connections

Each unit is equipped with a printer port and an RS-232C port. Simply connect the OCMA-350 to your printer or PC for remote control of the unit and ready access to measurement data.



The OCMA-350 has various special features to assure accurate measurements.

Stable readout

The unit's data evaluation function analyzes the readings, determines when the value has stabilized, and displays that stable value.

Self-monitoring

Messages on the LCD screen inform you of any electrical malfunctions, irregularities in measurement, and deterioration of parts, so that you can always be assured of proper functioning.

Simplified recordkeeping


The OCMA-350 keeps a record of the time and date of measurement along with each data set it records. When printing out, it provides you with both the time and the measurement, allowing for easy recordkeeping.




Measurement of residual oil

The OCMA-350 makes it easy to measure residual oil on pre-cleaned parts.


1 Choose a reasonably sized sample



2 Measure out an appropriate amount of solvent, and immerse the parts in the solvent until all the residual oil has dissolved.



3 Introduce the resulting solution to the measuring cell, and set the cell in place.




4 Press the start button and wait for the results.

HORIBA's S-316 Solvent

The OCMA-350 uses the highly effective, environmentally-safe S-316 extraction solvent to extract the oil components from oily water samples, soil samples, or product surfaces. The extract is analyzed using IR absorbance, a non-dispersive infrared spectrophotometric technique which is specific to hydrocarbons such as oil.


The OCMA-350 measures absorption in the 3.4 to 3.5 micrometer range. The two graphs at the right show the absorption spectra of (1) petroleum and (2) HORIBA's S-316 solvent. All hydrocarbons, including oils, absorb infrared radiation at about 3.4 to 3.5 micrometers. As a result, the unit can measure any hydrocarbons in the extraction solvent quickly and accurately, with no distortion of values due to the presence of the solvent.

HORIBA's S-316 Solvent can be recycled with the aid of the optional SR-300 Solvent Reclaimer. Recycling the solvent not only cuts your solvent costs, but also helps protect the environment.



Dimensional Outlines

Unit: mm (in)



Specifications

Application: Oil contamination of fresh and salt water, soil analysis

Principle: Solvent extraction, NDIR analysis (infrared spectroscopy)

Detector: Pyroelectric sensor

Cell: Quartz cell, cell length 20 mm

Measuring range and unit: 0 to 200 mg/l, 0 to 1000 mg/kg, 0 to 1 Abs.

Resolution:

mg/l: 0 to 99.9 mg/l, 0.1 mg/l
100 to 200 mg/l, 1 mg/l

mg/kg: 0 to 9.99 mg/kg, 0.01 mg/kg
10.0 to 99.9 mg/kg, 0.1 mg/kg
100 to 1000 mg/kg, 1 mg/kg

Abs: 0 to 1000 Abs, 0.001 Abs

Repeatability:

mg/l: 0 to 9.9 mg/l, ± 0.4 mg/l ± 1 digit
10.0 to 99.9 mg/l, ± 2.0 mg/l ± 1 digit
100 to 200 mg/l, ± 4 mg/l ± 1 digit

Abs: $\pm 1\%$ F.S.

Measurement: Manually controlled

Calibration: Automatic calibration after the calibration standard is introduced to the instrument

Extraction solvent: HORIBA S-316 solvent
Sample/solvent volume: Approx. 6.5 ml, ratio of sample to solvent is 1:1

Display: Measured value; 3 1/2 digits LCD with back-light
Message; Character display LCD with back-light (16 x 2 characters)

Functions: Self diagnostics, Auto hold function, Calendar clock

Output: RS-232C and centronics printer port

Ambient temperature: 0 to 40°C

Power requirement: 100 to 240 V AC $\pm 10\%$, 50/60 Hz, 60 VA

Dimensions: 200 (H) x 250 (W) x 285 (D) mm
7.9 (H) x 9.8 (W) x 11.0 (D) in

Weight: Approx. 5 kg/11 lb

Model SR-300 Solvent Reclaimer (optional)

Highly efficient, the SR-300 solvent reclaimer can reduce your per test solvent cost by up to 90%. This optional unit, designed especially for reclaiming HORIBA's S-316 solvent, features a double column of activated carbon and activated aluminum. It has a large filtering capacity, is easy to operate, and requires no electricity.

External dimensions:

200 (W) X 200 (D) X 600 (H) mm
7.9 (W) X 7.9 (D) X 23.6 (H) in



Horiba continues contributing to the preservation of the global environment through analysis and measuring technology.



Please read the operation manual before using this product to assure safe and proper handling of the product.

For more information, please contact us:

ExpotechUSA

10700 Rockley Road
Houston, Texas 77099
USA

281-496-0900 [voice]

281-496-0400 [fax]

E-mail: sales@expotechusa.com

Website: www.ExpotechUSA.com