


# Conductivity Meters

## Conductivity

Conductivity is a parameter used to measure the ionic concentration and activity of a solution. The more salt, acid or alkali in a solution, the greater its conductivity. The unit of conductivity is  $S/m$ , often also  $S/cm$ . The scale for aqueous solutions begins with pure water at a conductivity of  $0.05 \mu S/cm$  ( $77^\circ F/25^\circ C$ ). Naturally occurring waters such as drinking water or surface water have a conductivity in the range  $100 - 1000 \mu S/cm$ . At the upper end of the chart some acids and alkalines can be found.

Conductivity measurements are used for applications such as in the production of ultrapure water or determining the salinity of sea water.

Conductivity is measured by making a measurement of the electrical resistance. The simplest kind of measuring cell used consists of two similar electrodes. An alternating voltage applied to one of the electrodes causes the ions in the solution to migrate towards the electrodes. The more ions in the solution, the greater the current which flows between the electrodes. The instrument measures the current and uses Ohm's law to calculate first the conductance of the solution and then – by taken the cell data into account – the conductivity.



● recommended by WTW ○ conditionally applicable – not recommended

Application Range	inoLab®			Profi-Line	VARIO	Handheld Meters		
	Cond 720	Cond 730	Cond 740	Cond 197i	C <sub>cond</sub>	Cond 315i	Cond 330i	Cond 340i
Routine measurement	●	-	-	-	●	●	●	-
Routine measurement with documentation	-	●	●	●	-	-	-	●
AQA with documentation	-	●	●	●	-	-	-	●
R&D high precision	-	●	●	●	-	-	●	●
Control measurements	-	●	●	●	●	-	●	●
LIMS connection	-	●	●	●	-	-	-	○
Quality assurance	-	●	●	●	-	-	●	●
Training	●	●	●	○	●	●	●	○
Service	-	-	-	●	●	●	●	●
Laboratory measurements	●	●	●	●	●	-	-	○
Field measurements	-	-	-	●	-	●	●	●
Depth measurements	-	-	-	●	-	-	-	-
External control/PC connection/ PC control	-	●/●/-	●/●/●	●/●/-	-	-	-	●/●/-
Salinity/TDS measurement	●	●	●	●	●	only SAL	●	●
Specific resistance	●	●	●	-	-	●	●	●
Suitable for USP 27	●	●	●	●	-	-	●	●
Measurement of ultrapure water	●	●	●	●	●	●	●	●
Trace conductivity	●	●	●	●	-	-	●	●

see page

for conductivity measurements with multi-parameter instruments see page 50

Application Range Sensors	TetraCon® 325	TetraCon® 325/S	LR 325/01	LR325/001	TetraCon® DU/T	TA 197 LF
USP 27	-	-	●	●	-	-
Pharmaceutical water	○	-	●	●	-	-
Chemical water	○	-	-	-	●	-
Ground water	●	-	○	-	-	●
Surface waters	●	-	-	-	-	-
Depth measurements (barrages)	○	-	-	-	-	●
Laboratory measurements	●	-	●	●	-	-
Foods industry (juices)	●	-	-	-	○	-
Swimming pools	●	-	-	-	○	-
Pharmaceutics	●	-	●	○	○	-
Cosmetics/Detergents	○	●	-	-	-	-
Semi-conductor industry	-	-	●	●	-	-
Paint/Varnish (soluble)	●	○	-	-	-	-
Electroplating	●	-	-	-	-	-

applicable instruments:

all/except  
VARIO

all/except  
VARIO/  
315i

all/except  
VARIO

all/except  
VARIO/  
315i

all/except  
VARIO/  
315i

Cond  
197i



Parameter

pH

ORP

ISE

Oxygen  
(D.O.)

Conductivity

Multi-  
parameter

BOD/  
Respiration

Photometers

Turbidity

Colony  
Counter

Software/  
Printers

For more information, please contact us:

ExpotechUSA

10700 Rockley Road  
Houston, Texas 77099  
USA

281-496-0900 [voice]

281-496-0400 [fax]

E-mail: sales@expotechusa.com

Website: www.ExpotechUSA.com