

# ORP measurements

## ORP measurements

Reduction and oxidation are two central chemical terms that describe the ability of chemical agents to accept (reduction) or donate electrons (oxidation). In aqueous solutions, the ORP (Oxidation-Reduction Potential) voltage can be measured using a standard hydrogen electrode as reference. The reducing or oxidizing properties of a solution first are a matter of the reactants. By using an ORP electrode this change in potentials would be recorded as a positive or negative voltage.

ORP measurements monitor chemical reactions such as checking the denitrification of wastewater and disinfectant effect of detergents or the strength of plating baths. Measurement of ORP voltage is carried out with ORP combination electrodes. Similar to pH electrodes, these consist of a measuring electrode and a reference electrode. A metal electrode (usually made of precious metal such as platinum) is used in ORP combination electrodes in place of a glass membrane for carrying out the measuring function. The tendency for the chemical agents to accept or donate electrons determines the potential of the metal and thus the electrical potential of the combination electrode. ORP combination electrodes in use today contain a silver/silver chloride reference electrode, i. e. the indicated potential refers to this potential. Conversion to the standard hydrogen electrode system ( $U_H$ ) and that of the silver/silver chloride reference electrode is easily possible.

**SenTix® ORP reference electrode potential against the standard hydrogen electrode**

Temperature in °F (°C)	Potential in mV
32 (0)	+ 224
41 (5)	+ 221
50 (10)	+ 217
59 (15)	+ 214
68 (20)	+ 210
77 (25)	+ 207
86 (30)	+ 203
95 (35)	+ 200
104 (40)	+ 196
113 (45)	+ 192
122 (50)	+ 188
131 (55)	+ 184
140 (60)	+ 180
149 (65)	+ 176
158 (70)	+ 172

$$U_H = U_{Meas} + U_{Ref}$$

ORP measurements can be carried out with all WTW pH/mV meters.

## ORP measurements



Accurate ORP measurements can be made with any WTW pH/mV meters.



### Techn. Data and Ordering Information

	SenTix® ORP 103 648
Working temperature range	32 ... 212 °F (0 ... 100 °C)
Reference electrolyte	KCl 3 mol/l, Ag <sup>+</sup> -free
Sensor	Platinum
Sensor shape	Round, 0.16 in. (4 mm) dia.
Diaphragm	Ceramic
Shaft material	Glass
Shaft length / Shaft diameter	4.72 in. ±0.04 (120 mm ±2)/0.47 in. ±0.02 (12 mm ±0,5)
Temperature sensor	–
Connection	Plug head
Electrode cable (not included)	AS/DIN, AS/DIN-3 or AS/BNC
Electrode plug	DIN plug or BNC, as selected

### Ordering Informations for accessories

Testing and maintenance supplies for ORP measurements		Order No.
SORT/RH	Reagents for regenerating ORP electrodes consisting of activation powder (10 g) and clorina powder (30 g)	109 730
RH 28	ORP buffer solution 1 bottle of 250 ml: pH 7, U <sub>H</sub> = 427 mV	109 740

Parameter

pH

ORP

ISE

Oxygen  
(D.O.)

Conductivity

Multi-  
parameter

BOD/  
Respiration

Photometers

Turbidity

Colony  
Counter

Software/  
Printers

For more information, please contact us:

ExpotechUSA  
10700 Rockley Road  
Houston, Texas 77099  
USA

281-496-0900 [voice]

281-496-0400 [fax]

E-mail: sales@expotechusa.com

Website: www.ExpotechUSA.com