

60-MINUTE TIMER INSTRUCTIONS

SPECIFICATIONS

Display: 2¹/₂"-diameter dial
Range: 0 to 60 minutes
Resolution: 1 minute
Alarm: Audible
Size: 3¹/₄H x 3W x 2D
Weight: 4.8 oz

OPERATION

Turn the knob clockwise to 60 and then turned back to the desired time.

The bell will ring at the end of the selected time.

For more information, please contact us:

[ExpotechUSA](#)
[10700 Rockley Road](#)
[Houston, Texas 77099](#)
[USA](#)

[281-496-0900 \[voice\]](#)

[281-496-0400 \[fax\]](#)

E-mail: sales@expotechusa.com

Website: www.ExpotechUSA.com