



## Performance Specifications of the 3100 Canister Cleaner

- Cleans up to 16 6L canisters or 63 MiniCans at a time.
- Cycles between evacuation and fill with humid zero air or N(g).
- Oilless pumps eliminate the need for cryogenic or sorbent oil vapor traps.
- Choose band heaters or a 6Lx6 position oven.
- All stainless tubing is internally electropolished.
- Optional Silonite fused silica lining of manifold and vacuum system virtually eliminates surface absorption while decreasing canister cleanup times.
- Cleans and leak checks canisters within 1-3 hours.
- Incorporates separate high vacuum and rough pumps for improved efficiency.
- Achieve < 20 mtorr vacuum on final pump-down with Entech supplied vacuum system.
- Automated leak checking capabilities and QAQC reporting
- Controlled using Smartlab software operating on a Windows (3.1, 95) PC.
- Fill/evacuation cycles controlled by time or pressure setpoints.
- Compatible with Entech Silonite Coated Canisters, SUMMA Canisters, and other Spherical Cans.
- Compatible with Entech MiniCans with Miniature Quick Connect for lower cost and ease of use.



For more information on the above products and services:

Please contact us:

ExpoTechUSA  
10700 Rockley Road  
Houston, Texas 77099  
USA

281-496-0900 [voice]

281-496-0400 [fax]

E-mail: [sales@expotechusa.com](mailto:sales@expotechusa.com)

Web: [www.ExpoTechUSA.com](http://www.ExpoTechUSA.com)