

Laboratory Safety Ventilation

*Where You Need It,
When You Need It*



The *Labconco Express* program speeds delivery of selected safety ventilation products to you. After receipt of the order, products are shipped in 48 hours or less.

We've selected our most popular safety ventilation products for the *Express* program. The Protector® Fiberglass Hood heads the list, with 4', 5' and 6' models available for rapid delivery. We've included four factory installed service fixtures and duplex electrical receptacles on each hood to eliminate the time and expense involved with field installation of these accessories — another time-saving feature.

Where space is at a premium, the Basic™ 47 Hood is an economical alternative to meet your ventilation requirements. Its 53" height (including built-in blower) and 25" depth make it a perfect fit for trailers and temporary buildings.

Hoods are available with or without built-in blowers. Coated Steel Blowers for hoods without them are also offered. Your choice of three different base cabinet styles and companion work surfaces complete your installation.

Select the Paramount® Filtered Enclosure for protection *without ducting* from organic vapors, acid gases, formaldehyde odors or ammonia fumes. The Paramount houses Color-Smart™ Filter Cells which adsorb or treat

fumes and vapors returning filtered air to the laboratory. Select the Paramount Duet™ Filtered Enclosure, which used both carbon and HEPA filter cells, for applications involving gaseous and particulate contaminants. Filter Cells are sold separately. A Mobile Base Stand with Dished Work Surface is also available for rapid delivery.

The Purifier® Delta™ Series Safety Cabinet provides protection for both the product and the operator from particulate contamination including blood-borne pathogens and antineoplastic drugs. The 4' Class II Type A/B3 cabinet includes a service fixture, ultraviolet light and two electrical duplexes with splash covers inside the work area. An accessory stand is adjustable to various heights.

For personnel protection from particulates that do not require product protection, the Purifier® HEPA Filtered Enclosure is ideal. Applications include weighing powders, asbestos handling and media preparation. This enclosure is offered in 2' and 3' widths for placement atop an accessory Dished Work Surface and existing cabinetry.

To order, or for additional information on any *Labconco Express* product, contact your laboratory supply dealer or Labconco Customer Service at 800-821-5525 or 816-333-8811.

Express Products



Basic 47 Hood #22473-00 with Work Surface #48828-03 and Acid Storage Cabinet #99010-00

Basic™ 47 Laboratory Hood, Base Cabinets and Flat Work Surface

- 22473-00 Basic 47 Hood with blower
- 99010-00 4' Acid Storage Cabinet (*requires work surface*)
- 99020-00 4' Solvent Storage Cabinet (*requires work surface*)
- 99000-00 4' Standard Base Cabinet (*requires work surface*)
- 48828-03 4' Flat Solid Epoxy Work Surface

Protector® 48 Fiberglass Hoods, Blower, Base Cabinets and Dished Work Surface

- 48803-00 Protector 48 Fiberglass Hood with blower and four fixtures
- 48804-00 Protector 48 Fiberglass Hood without blower, with four fixtures (*requires blower*)
- 70682-00 Coated Steel Blower for roof mounting
- 99010-00 4' Protector Acid Storage Cabinet (*requires work surface*)
- 99020-00 4' Protector Solvent Storage Cabinet (*requires work surface*)
- 99000-00 4' Standard Base Cabinet (*requires work surface*)
- 48835-00 4' SpillStopper Solid Epoxy Dished Work Surface



Express Products

Protector® 60 Fiberglass Hoods, Base Cabinets and Dished Work Surface

- 60803-00 Protector 60 Fiberglass Hood with blower and four fixtures
 - 60804-00 Protector 60 Fiberglass Hood without blower, with four fixtures (requires blower)
 - 70688-00 Coated Steel Blower for roof mounting
 - 99012-00 30" Protector Acid Storage Cabinet (requires work surface)
 - 99022-00 30" Protector Solvent Storage Cabinet (requires work surface)
 - 99002-00 30" Standard Base Cabinet (requires work surface)
- Any combination of two base cabinets listed above may be ordered.*
- 60835-00 5' SpillStopper Solid Epoxy Dished Work Surface

Protector® 72 Fiberglass Hoods, Base Cabinets and Dished Work Surface

- 72803-00 Protector 72 Fiberglass Hood with blower and four fixtures
 - 72804-00 Protector 72 Fiberglass Hood without blower, with four fixtures (requires blower)
 - 70688-00 Coated Steel Blower for roof mounting
 - 99011-00 36" Protector Acid Storage Cabinet (requires work surface)
 - 99021-00 36" Solvent Storage Cabinet (requires work surface)
 - 99001-00 36" Standard Base Cabinet (requires work surface)
- Any combination of two base cabinets listed above may be ordered.*
- 72835-00 6' SpillStopper Solid Epoxy Dished Work Surface

Paramount® Filtered Enclosures, Filter Cells and Mobile Base Stand

- 69101-00 3' Paramount Filtered Enclosure (requires one package of Carbon-Based Filter Cells)
- 69101-10 3' Paramount Duet™ Filtered Enclosure (requires one package of Carbon-Based Filter Cells and one package of HEPA Filter Cells)
- 69102-00 Color-Smart Ammonia Carbon-Based Filter Cells, package of 2
- 69103-00 Color-Smart Acid Gas Carbon-Based Filter Cells, package of 2
- 69104-00 Color-Smart Organic Carbon-Based Filter Cells, package of 2
- 69105-00 Color-Smart Formaldehyde Carbon-Based Filter Cells, package of 2
- 69332-00 Color-Smart HEPA Filter Cells, package of 2
- 69107-00 3' Mobile Base Stand with solid epoxy dished work surface

4' Purifier® Delta™ Series Class II, Type A/B3 Safety Cabinet and Base Stand

- 36209-04 4' Purifier Delta Series Safety Cabinet with UV light and one service fixture
- 37304-00 4' Base Stand, adjustable from 27 1/2- 34 1/2" high

Purifier HEPA Filtered Enclosures and Work Surfaces

- 37200-00 2' Purifier HEPA Filtered Enclosure
- 37300-00 3' Purifier HEPA Filtered Enclosure
- 48621-00 2' Dished Epoxy Work Surface
- 48631-00 3' Dished Epoxy Work Surface



Protector 48 Hood #48803-00 with Work Surface #48835-00 and Acid Storage Cabinet #99010-00



Purifier Delta Series Class II Safety Cabinet #36209-04 with Base Stand #37304-00



Paramount Filtered Enclosure #69101-00



Purifier HEPA Filtered Enclosure #37300-00 with Dished Epoxy Work Surface #48631-00



For more information, please contact us:

[ExpotechUSA](#)
[10700 Rockley Road](#)
[Houston, Texas 77099](#)
[USA](#)

[281-496-0900 \[voice\]](#)

[281-496-0400 \[fax\]](#)

E-mail: sales@expotechusa.com

Website: www.ExpotechUSA.com