



# Laboratory Carts & Benches



*Protecting your  
laboratory environment*

**LABCONCO**

Labconco Laboratory Carts and Benches fit in every lab. These versatile carts have endless uses from transporting equipment and supplies to providing auxiliary bench space. Labconco Carts and Benches feature:

- Ready-to-roll, full factory assembly.
- Toe-locking casters.
- Ergonomic, extruded rubber handle grips.
- Durable glacier white, epoxy-coated tubular steel frames.
- ID tag to ensure your cart makes a round trip back to you.
- Lifetime guarantee – if your Labconco Cart or Bench does not live up to your expectations, we will replace it FREE.



## Flexi-Bin™ Carts

- For safe transport of glassware and chemical containers.
- Removable, corrosion-resistant polyethylene bin, 22.25 x 17.00" (56.5 x 43.2 cm), is included. Inside dimensions: 20.25 x 15.25" (51.4 x 38.7 cm). Bins are offered in two overall depths of 4" and 10" with inside depths of 3.8" and 9.8", and are interchangeable.
- Include three long (20.25") and four short (15.25") polyethylene dividers for partitioning the bin to securely hold items in upright positions.
- Bottom open shelf holds additional supplies.
- 4" diameter casters, two with toe locks.
- Lifetime guarantee.
- Dimensions: 34.88 x 19.00 x 36.38" high (88.6 x 48.3 x 92.4 cm)

Catalog #	Description	Shipping Weight
8010000	Flexi-Bin Cart with 4" bin	44 lbs. (20 kg)
8010001	Flexi-Bin Cart with 10" bin	44 lbs. (20 kg)
8068000	Bin, 4" deep, with dividers	5 lbs. (2.3 kg)
8068100	Bin, 10" deep, with dividers	10 lbs. (4.5 kg)

## Stockroom Cart



#8021000

- Catalog #8021000.
- Step ladder and cart in one.
- Two removable, watertight, epoxy-coated steel pans with gray vinyl mats, and one vinyl-coated wire basket.
- Sturdy, two-step ladder with non-skid surface conforms to OSHA and ANSI design guidelines.
- Spring-loaded, oversized steps with skid-resistant, rubber feet prevent cart from moving while user stands on the ladder.
- Steps withstand loads up to 1000 lbs. (453.6 kg).
- 4" diameter casters, two with toe locks.
- Lifetime guarantee.
- Dimensions: Cart is 19.00 x 39.50 x 43.00" high (48.3 x 100.3 x 109.2 cm). Top pan is 6.50 x 16.75 x 2.00" (16.5 x 42.5 x 5.1 cm). Bottom pan is 17.50 x 16.75 x 2.00" (44.5 x 42.5 x 5.1 cm). Basket is 16.00 x 18.00 x 4.00" (40.6 x 45.7 x 10.2 cm).
- Shipping weight 69 lbs. (31.3 kg).

## Utility Cart

- Catalog #8007000.
- Two open shelves transport and support small electrical equipment. May be used as a compact work station for laboratory equipment.
- Power strip with 115 volts, 60 Hz, 12 amps capacity, has 6 electrical outlets, cord bracket and a 3-wire cord with plug.
- Two corrosion-resistant, one-piece, vacuum-formed ABS plastic drawers, 15.75 x 14.88 x 4.20" high (40.0 x 37.8 x 10.7 cm), hold small labware and supplies.
- 4" diameter casters, two with toe locks.
- Lifetime guarantee.
- Dimensions: 34.88 x 19.00 x 36.38" high (88.6 x 48.3 x 92.4 cm)
- Shipping weight 44 lbs. (20kg)



#8007000

The 18-gauge welded steel frames of 1" square steel tubing are painted with tough, corrosion-resistant, dry powder epoxy.

Benches have 2" square steel tubing.

# Durable

# Quiet

Four-inch diameter casters absorb vibration over rough, uneven floors. Two toe-locks keep carts stationary when needed. Benches have 5" diameter, non-marking toe-locking casters with bearings for quiet operation.

## Portable Table

- Catalog #8025000.
- Labconco's most popular cart has two open shelves that provide easy transfer of equipment from countertop to cart.
- 4" diameter casters, two with toe locks.
- Lifetime guarantee.
- Dimensions: 34.88 x 19.00 x 36.38" high (88.6 x 48.3 x 92.4 cm).
- Shipping weight 44 lbs. (20 kg).



## Chemical Cart

- Catalog #8020000.
- Similar to the Portable Table but rigid PVC extrusions on two sides prevent containers from sliding off. The lower shelf also has rigid plastic guard rails on both ends.
- 4" diameter casters, two with toe locks.
- Lifetime guarantee.
- Dimensions: 34.88 x 19.00 x 36.38" high (88.6 x 48.3 x 92.4 cm).
- Shipping weight 46 lbs. (20.9 kg).



#8045000



#8032500

## Glassware Carts

- Carry glassware, plasticware, even small instruments, in vinyl-coated wire baskets.
- Choose 2-basket or 4-basket style.
- Complete with removable plastic drip pan to catch spills.
- 4" diameter casters, two with toe locks.
- Lifetime guarantee.
- Dimensions: 34.88 x 19.00 x 36.38" high (88.6 x 48.3 x 92.4cm).

Catalog #	Description	Shipping Weight
8032500	2-Basket Glassware Cart. Each basket measures 32.00 x 18.00 x 4.00" (81.3 x 45.7 x 10.2 cm)	43 lbs. (19.5 kg)
8045000	4-Basket Glassware Cart. Each basket measures 16.00 x 18.00 x 4.00" (40.6 x 45.7 x 10.2 cm)	43 lbs. (19.5 kg)
8040100	Large Basket, 32.00 x 18.00 x 4.00" (81.3 x 45.7 x 10.2 cm)	5 lbs. (2.3 kg)
8040200	Small Basket, 16.00 x 18.00 x 4.00" (40.6 x 45.7 x 10.2 cm)	2 lbs. (0.9 kg)

## Stainless Steel Cart

- Catalog #8019000.
- Holds a variety of laboratory supplies, equipment and miscellaneous containers.
- Two removable polished Type 304 stainless steel pans are 2" (51 mm) deep on three sides with one open side for easy equipment transfer from countertop to cart.
- 4" diameter casters, two with toe locks.
- Lifetime guarantee.
- Dimensions: 34.88 x 19.13 x 36.38" high (88.6 x 48.6 x 92.4 cm).
- Shipping weight 56 lbs. (25.4 kg).



#8019000



#8052500

## Instrument Desk

- Catalog #8052500.
- Serves as a desk, instrument table or laboratory cart.
- Support bar is centered to permit comfortable seating from either side.
- Built to standard desk height.
- Work surface is phenolic board with a melamine high-pressure laminate surface that resists stains and corrosion.
- 4" diameter casters, two with toe locks.
- Lifetime guarantee.
- Dimensions: 42.13 x 19.00 x 31.63" high. (107 x 48.3 x 80.3 cm).
- Shipping weight 35 lbs. (15.9 kg).

## Pan Cart

- Catalog #8047500.
- Two 2" (51 mm) deep epoxy-coated steel, removable watertight pans transport wet or soiled glassware and other labware.
- Removable gray rubber mats line pans, prevent loads from shifting.
- 4" diameter casters, two with toe locks.
- Lifetime guarantee.
- Dimensions: 34.88 x 19.00 x 36.38" high (88.6 x 48.3 x 92.4 cm).
- Shipping weight 58 lbs. (26.3 kg).



#8047500



#8000300

## Carboy Caddy

- Catalog #8000300.
- Safe way to move large carboys or cartons.
- Helps eliminate heavy lifting and back strain.
- 3" diameter casters provide low profile, elevating shelf only 4.25" above the floor.
- Shelf is made of phenolic board with a melamine high-pressure laminate surface.
- Handle measures 38.50" high and 19.00" wide (97.8 x 48.3 cm).
- Shelf dimensions: 20.00 x 23.63" (50.8 x 60 cm).
- Supports loads up to 360 lbs. (163.3 kg).
- Lifetime guarantee.
- Shipping weight 21 lbs. (9.5 kg).

# Strong

Labconco Carts support heavy loads.  
Carts hold up to 400 pounds; Benches up to 540.



#8060000

## Mobile Bench

- Catalog #8060000.
- Heavy duty auxiliary bench may be transported from lab to lab.
- 5" diameter, non-marking polyurethane casters with toe locks.
- Lifetime guarantee.
- Dimensions: 48 x 24 x 37 3/16" high. (121.9 x 61 x 94.5 cm). Work surface is 1.13" (29 mm) thick.
- Shipping weight 111 lbs. (50.3 kg).

## Variable Height Bench

- Catalog #8075000.
- Heavy duty auxiliary bench has adjustable height.
- 5" diameter, non-marking polyurethane casters with toe locks.
- Four adjustments from 29.63 to 37.38" (75.1 to 94.9 cm) using steel pins with nuts inserted in vertical frame members.
- Lifetime guarantee.
- Work surface dimensions: 38.00 x 28.00 x 1.13" thick (96.5 x 71.1 cm x 29 mm).
- Shipping weight 111 lbs. (50.3 kg).



#8075000

# Flexible

Hydraulic Lift Base Stands may be easily raised and lowered with the turn of a crank or push of a button.



#3780400 shown with #4882805 and #3784000.

**3784000 Caster Kit.** Includes four 5" diameter, non-marking polyurethane casters with toe locks. Load rating for each caster is 250 pounds. Casters have bearings for quiet operation. Casters replace the fixed feet that are included with the stand. **Installation is required.** Shipping weight 10 lbs. (4.5 kg)

\* Americans with Disabilities Act  
\*\* International electrical configuration

## Hydraulic Lift Base Stands

- When coupled with a work surface, the stands may be used as an adjustable height table. Work Surfaces are sold separately.
- Durable 1.75" tubular, epoxy-coated steel frame with rear crossbar supports. Sub-frame is 2.0" x 3.0" .
- Working height is infinitely adjustable from 26.75" to 34.75" high.
- Four powder-coated, die-cast feet with levelers.
- Supports loads up to 1000 pounds.
- Optional accessory 5" diameter, toe-locking Caster Kit is available for field installation.
- ADA\*-compliant.

### Flat Solid Epoxy Work Surfaces

Nominal Width	Catalog #	Overall Width x Depth x Thickness	Shipping Weight
3 feet	4882805	36.0" x 30.0" x 1.25"	90 lbs. (41 kg)
4 feet	4882803	48.0" x 30.0" x 1.25"	110 lbs. (50 kg)
6 feet	4882804	72.0" x 30.0" x 1.25"	205 lbs. (93 kg)

### Manual Hydraulic Lift Base Stands

Nominal Width	Catalog #	Overall Width x Depth	Shipping Weight
3 feet	3780400	36.0" x 29.0"	100 lbs. (45 kg)
4 feet	3780401	48.0" x 29.0"	120 lbs. (54 kg)
6 feet	3780402	72.0" x 29.0"	150 lbs. (68 kg)

### Electric Hydraulic Lift Base Stands

Nominal Width	Catalog #	Overall Width x Depth	Electrical Requirements	Shipping Weight
3 feet	3780300	36.0" x 29.0"	115 volts, 60 Hz, 2.4 amps	100 lbs. (45 kg)
3 feet	3780303**	36.0" x 29.0"	230 volts, 50 Hz, 0.8 amp	100 lbs. (45 kg)
4 feet	3780301	48.0" x 29.0"	115 volts, 60 Hz, 2.4 amps	120 lbs. (54 kg)
4 feet	3780304**	48.0" x 29.0"	230 volts, 50 Hz, 0.8 amp	120 lbs. (54 kg)
6 feet	3780302	72.0" x 29.0"	115 volts, 60 Hz, 2.4 amps	150 lbs. (68 kg)
6 feet	3780305**	72.0" x 29.0"	230 volts, 50 Hz, 0.8 amp	150 lbs. (68 kg)



Protecting your  
laboratory environment

**LABCONCO**

For more information, please contact us:

ExpotechUSA

10700 Rockley Road  
Houston, Texas 77099  
USA

281-496-0900 [voice]

281-496-0400 [fax]

E-mail: sales@expotechusa.com

Website: www.ExpotechUSA.com