



Class II, Type B2 Remote Blower - 3663500




Blower RPM	Motor Sheave Setting	CFM	Static Pressure	Motor Current (A)
2074	max	1179	3.98	6.0
2102	max	826	4.14	5.0
2123	max	518	4.38	4.3
1978	1 turn	1337	3.64	6.0
2008	1 turn	945	3.75	5.0
2036	1 turn	555	3.96	4.3
1892	2 turns	1448	3.34	6.0
1914	2 turns	1075	3.37	5.0
1938	2 turns	726	3.51	4.3
1812	3 turns	1565	3.08	6.0
1826	3 turns	1175	3.12	5.0
1854	3 turns	650	3.21	4.0
1705	4 turns	1736	2.50	6.0
1715	4 turns	1600	2.66	5.5
1728	4 turns	1311	2.79	4.9
1760	4 turns	959	2.93	4.4
1765	4 turns	736	2.97	4.0
1622	5 turns	1905	1.89	6.0
1632	5 turns	1729	2.18	5.6
1645	5 turns	1458	2.50	5.0
1665	5 turns	1138	2.64	4.5
1682	5 turns	807	2.71	4.0
1685	5 turns	618	2.74	3.8

Class II, Type A2 Remote Blower – 3668000
For 3' and 4' Enclosures



Blower RPM	Motor Sheave Setting	CFM	Static Pressure	Motor Current (A)
1188	max	576	0.15	4.2
1190	max	555	0.25	4.1
1196	max	448	0.50	3.9
1204	max	284	0.72	3.7
1208	max	163	0.84	3.7
1128	1 turn	555	0.14	4.0
1128	1 turn	530	0.25	3.9
1132	1 turn	453	0.42	3.8
1136	1 turn	312	0.59	3.7
1142	1 turn	146	0.74	3.5
1064	2 turns	530	0.13	3.9
1066	2 turns	490	0.25	3.8
1069	2 turns	407	0.40	3.7
1071	2 turns	292	0.53	3.6
1076	2 turns	136	0.66	3.5
1000	3 turns	497	0.11	3.8
1003	3 turns	443	0.26	3.7
1006	3 turns	314	0.42	3.6
1010	3 turns	134	0.58	3.6
937	4 turns	448	0.10	3.7
939	4 turns	395	0.24	3.6
941	4 turns	278	0.38	3.5
943	4 turns	119	0.50	3.5

**Class II, Type A2 Remote Blower – 3668001
For 6' Enclosures**



Blower RPM	Motor Sheave Setting	CFM	Static Pressure	Motor Current (A)
1503	max	760	0.25	8.3
1507	max	700	0.52	8.1
1514	max	609	0.77	8.0
1522	max	454	1.03	7.8
1534	max	193	1.35	7.6
1434	1 turn	717	0.23	8.1
1440	1 turn	661	0.48	8.1
1443	1 turn	576	0.72	7.9
1451	1 turn	418	0.97	7.8
1458	1 turn	185	1.22	7.7
1367	2 turns	688	0.21	8.1
1370	2 turns	638	0.43	8.0
1373	2 turns	550	0.65	7.9
1378	2 turns	397	0.87	7.8
1384	2 turns	178	1.10	7.7
1295	3 turns	649	0.18	7.9
1299	3 turns	567	0.46	7.8
1303	3 turns	423	0.71	7.7
1309	3 turns	174	0.98	7.6
1224	4 turns	605	0.17	7.8
1226	4 turns	536	0.41	7.7
1230	4 turns	395	0.63	7.6
1236	4 turns	160	0.87	7.6

For more information, please contact us:

[ExpotechUSA](#)
[10700 Rockley Road](#)
[Houston, Texas 77099](#)
[USA](#)

[281-496-0900 \[voice\]](#)

[281-496-0400 \[fax\]](#)

E-mail: sales@expotechusa.com

Website: www.ExpotechUSA.com