

Purifier Replacement Filter Spec Sheet

Cabinet Model #	Cabinet Description	Filter Type	Dimensions (")	P/N	Diff Press (" wq)	Rated CFM
36204/36205	3' Class II Type A2 Purifier Delta	Supply HEPA	18 x 36 x 3.06	3739400	0.67	335
		Exhaust HEPA	18 x 12 x 11.5	3739500	0.51	266
36208/36209	4' Class II Type A2 Purifier Delta	Supply HEPA	18 x 48 x 3.06	3739401	0.67	445
		Exhaust HEPA	18 x 24 x 11.5	3739501	0.51	354
36212/36213	6' Class II Type A2 Purifier Delta	Supply HEPA	18 x 72 x 3.06	3739402	0.67	665
		Exhaust HEPA	18 x 30 x 11.5	3739502	0.51	529
36210	4' Class II Type B2 Purifier Delta	Supply HEPA	18 x 48 x 3.06	3739401	0.67	445
		Exhaust HEPA	18 x 30 x 5.88	3759500	0.69/0.79	770
36214	6' Class II Type B2 Purifier Delta	Supply HEPA	18 x 72 x 3.06	3739402	0.67	665
		Exhaust HEPA	18 x 48 x 5.88	3759600	0.64/0.74	1175
37200-	2' Class I/HEPA Filtered Enclosure	Exhaust HEPA	18 x 18 x 3.06	3707900	0.60	189
37300-	3' Class I/HEPA Filtered Enclosure	Exhaust HEPA	18 x 30 x 3.06	3707901	0.60	327
37500-03/-04/-05	3' PVC Total Exhaust Clean Bench, Trace Metals Work Station, Forensic Enclosure, Purifier Class 100 Chemical Station	Exhaust HEPA	18 x 30 x 3.06	3725900	0.60	375
		Prefilter	4.62 x 4.62 x 1.0	3727500	N/A	N/A
3740002/3740022	2' Purifier PCR Enclosure	Supply HEPA	18 x 18 x 3.06	3707900	0.60	189
		Prefilter	15 x 20.25 x .88	3718000	N/A	N/A
3750002/3750022	3' Purifier PCR Enclosure	Supply HEPA	18 x 30 x 3.06	3707901	0.60	327
		Prefilter	15 x 32.25 x .88	3718001	N/A	N/A
3740000/-01/-20/-21	2' Vertical Clean Bench	Supply HEPA	18 x 18 x 3.06	3707900	0.60	189
		Prefilter	15 x 20.25 x .88	3718000	N/A	N/A
3750000/-01/-20/-21	3' Vertical Clean Bench	Supply HEPA	18 x 30 x 3.06	3707901	0.60	327
		Prefilter	15 x 32.25 x .88	3718001	N/A	N/A
36000-	3' Horizontal Clean Bench	Supply HEPA	30 x 36 x 5.88	1489000	1.00	1135
		Prefilter	17.75 x 17.75 x .75	3768900	N/A	N/A
36125-	4' Horizontal Clean Bench	Supply HEPA	30 x 48 x 5.88	1489002	1.00	1525
		Prefilter	17.75 x 23.75 x .75	3768901	N/A	N/A
36100-	6' Horizontal Clean Bench	Supply HEPA	30 x 72x 5.88	1489001	1.00	2300
		Prefilter	17.75 x 23.75 x .75	3768901	N/A	N/A

For more information, please contact us:

[ExpotechUSA](#)
[10700 Rockley Road](#)
[Houston, Texas 77099](#)
[USA](#)

[281-496-0900 \[voice\]](#)

[281-496-0400 \[fax\]](#)

E-mail: sales@expotechusa.com

Website: www.ExpotechUSA.com