



Protecting your  
laboratory environment

**LABCONCO®**



## Product Specifications For **UNIVERSAL™ BASE STANDS**

**80882-00, 80882-01, 80883-00, 80883-01, 80860-00, 80860-01,  
80884-00, 80884-01, 80870-00, 80870-01, 80894-00, 80894-01,  
80873-10, 80873-11, 80873-12, 80873-13,  
80880-00, 80880-01, 80890-00, 80890-01, 80895-00, 80895-01**

## **GENERAL DATA**

### **APPLICATION**

Universal Base Stands support a variety of enclosures: Purifier® Class II Safety Cabinets, Purifier® HEPA Filtered Enclosures with Work Surfaces, Purifier® Class I Safety Enclosures with Work Surfaces, Purifier® Vertical Clean Benches with Work Surfaces, Purifier® Trace Metals Work Stations, Purifier® Forensic Enclosures, Purifier® Class 100 Chemical Stations, Purifier® Horizontal Clean Benches, Protector® VS Ventilation Stations with Work Surfaces, Protector® Laboratory Hoods with Work Surfaces, Basic™ Laboratory Hoods with Work Surfaces, and Protector® Work Stations with Work Surfaces. See chart on back for specific catalog numbers.

### **CONSTRUCTION**

Universal Base Stands are constructed of durable charcoal epoxy-coated steel with removable leveling feet. Models 80870-00, 80870-01, 80873-10, and 80873-11 include a 2 1/2" diameter cutout with plug on the right side to allow pass-through of plumbing

lines. Base Stands are fully factory assembled except for the following: 80873-10, 80873-11, 80873-12, and 80873-13 which **require assembly**.

### **SITTING HEIGHT**

The following Universal Base Stands are designed at sitting height, 28 3/4" (73 cm): 80882-00, 80883-00, 80860-00, 80884-00, 80870-00, 80873-10, 80873-12, 80880-00, 80890-00, 80894-00 and 80895-00.

### **STANDING HEIGHT**

The following Universal Base Stands are designed at standing height, 34 3/4" (88.3 cm): 80882-01, 80883-01, 80860-01, 80884-01, 80894-01, 80870-01, 80873-11, 80873-13, 80880-01, 80890-01 and 80895-01.


### **NSF INTERNATIONAL APPROVED**

All Universal Base Stands designed to support Purifier® Class II Safety Cabinets are approved by the National Sanitation Foundation International.

# DIMENSIONAL DATA

CATALOG NUMBER	NOMINAL WIDTH	A	B	C	FOR USE WITH:	SHIP. WT. lbs./kg
80882-00	2 feet	24"	28 3/4"	—	2' Purifier Class I Safety Enclosures, PCR Enclosures, HEPA Filtered Enclosures and Vertical Clean Benches, Protector VS Ventilation Stations	60/27
80882-01	2 feet	24"	34 3/4"	—		66/30
80883-00	3 feet	36"	28 3/4"	—	3' Purifier Class I Safety Enclosures, PCR Enclosures, HEPA Filtered Enclosures and Vertical Clean Benches, Protector VS Ventilation Stations	86/39
80883-01	3 feet	36"	34 3/4"	—		93/42
80860-00	3 feet	38 3/8"	28 3/4"	—	3' Purifier Class II Safety Cabinets and Horizontal Clean Benches	86/39
80860-01	3 feet	38 3/8"	34 3/4"	—		93/42
80884-00	4 feet	48"	28 3/4"	—	4' Protector VS Ventilation Stations	93/42
80884-01	4 feet	48"	34 3/4"	—		117/53
80870-00 80873-10*	4 feet	50 3/8"	28 3/4"	17 3/8"	4' Purifier Class II Safety Cabinets and Horizontal Clean Benches, Protector Work Stations	93/42 117/53
80870-01 80873-11*	4 feet	50 3/8"	34 3/4"	23 3/8"		99/45 123/56
80894-00	4 feet	47"	28 3/4"	—	Protector 48 Laboratory Hoods, Protector 96 Laboratory Hoods (2 base stands required), Basic 47 Laboratory Hoods	93/42
80894-01	4 feet	47"	34 3/4"	—		117/53
80880-00	5 feet	59"	28 3/4"	—	Protector 60 Laboratory Hoods	126/57
80880-01	5 feet	59"	34 3/4"	—		132/60
80890-00 80873-12*	6 feet	74 3/8"	28 3/4"	—	6' Purifier Class II Safety Cabinets and Horizontal Clean Benches	130/59 137/62
80890-01 80873-13*	6 feet	74 3/8"	34 3/4"	—		136/62 143/65
80895-00	6 feet	70"	28 3/4"	—	Protector 72 Laboratory Hoods, Basic 70 Laboratory Hoods	130/59
80895-01	6 feet	70"	34 3/4"	—		137/62

\* Unassembled



For more information, please contact us:

[ExpotechUSA](#)  
[10700 Rockley Road](#)  
[Houston, Texas 77099](#)  
[USA](#)

[281-496-0900 \[voice\]](#)

[281-496-0400 \[fax\]](#)

E-mail: [sales@expotechusa.com](mailto:sales@expotechusa.com)

Website: [www.ExpotechUSA.com](http://www.ExpotechUSA.com)