

6 REASONS

Why Labconco Carts & Benches Are Your Best Laboratory Transportation Value

LABCONCO
LABORATORY
CARTS & BENCHES
LIFETIME
GUARANTEE



1
DURABLE

Welded frame of 1" (25 mm) square steel tubing provides durability. Benches have 2" (51 mm) square tubing. Baked-on epoxy powder coating is abrasion and corrosion resistant.

2
STRONG

Capable of transporting heavy loads. Most cart models support up to 400 lbs. (181.4 kg). Benches hold up to 540 lbs. (245 kg).

3
QUIET

Four-inch diameter hard rubber casters are easy to steer and absorb vibration for smooth, quiet trips over rough, uneven floors. Benches have 5" (127 mm) diameter toe-locking casters with protective rubber guards.

4
ASSEMBLED

Shipped from our factory ready to roll.

If your Labconco Cart or Bench does not live up to your expectations, we will

GUARANTEED

5

replace it **FREE!**

6
12 MODELS

Each cart is specifically designed for the unique demands of the laboratory. Most carts are less than 20" wide to fit through narrow corridors and doorways.

LABCONCO

LABORATORY CARTS



#80475-00

PAN CART

- Two 2" (51 mm) deep epoxy-coated steel, removable watertight pans.
- Removable black rubber mats line pans, prevent loads from shifting.
- Dimensions: 34 7/8 x 19 x 36 3/8" high (88.6 x 48.3 x 92.4 cm).
- Shipping weight 58 lbs. (26.3 kg).

STAINLESS STEEL CART

- Holds a variety of laboratory supplies, equipment and miscellaneous containers.
- Two removable polished Type 304 stainless steel pans are 2" (51 mm) deep on three sides with one open side for easy equipment transfer from countertop to cart.
- Dimensions: 34 7/8 x 19 1/8 x 36 3/8" high (88.6 x 48.6 x 92.4 cm).
- Shipping weight 56 lbs. (25.4 kg).



#80190-00

INSTRUMENT DESK

- Serves as a desk, instrument table or laboratory cart.
- Support bar is centered to permit comfortable seating from either side.
- Built to standard desk height.
- Toe locks on two casters secure cart for stationary use.
- Work surface is phenolic board with a melamine high-pressure laminate surface which resists stains and corrosion.
- Dimensions: 42 1/8 x 19 x 31 5/8" high. (107 x 48.3 x 80.3 cm).
- Shipping weight 35 lbs. (15.9 kg).



#80525-00



#80003-00

LABORATORY CADDY

- Safe way to move large carboys or cartons.
- Helps eliminate heavy lifting and back strain.
- Holds up to 360 lbs. (163.3 kg).
- 3" (76 mm) diameter casters provide low profile, elevating shelf only 4 1/4" (114 mm) above the floor.
- Shelf is made of phenolic board with a melamine high-pressure laminate surface.
- Handle measures 38 1/2" high and 19" wide (97.8 x 48.3 cm).
- Shelf dimensions: 20 x 23 5/8" (50.8 x 60 cm).
- Shipping weight 21 lbs. (9.5 kg).

CHEMICAL CART

- Similar to the Portable Table below but with guard rails to prevent containers from sliding off.
- Dimensions: 34 7/8 x 19 x 36 3/8" high (88.6 x 48.3 x 92.4 cm).
- Shipping weight 46 lbs. (20.9 kg).



#80450-00 4-Basket Cart



GLASSWARE CARTS

- Carry glassware, plasticware, even small instruments, in vinyl coated wire baskets.
- Choose 2-basket or 4-basket style.
- Complete with removable plastic drip pan to catch spills.
- Dimensions: 34 7/8 x 19 x 36 3/8" high (88.6 x 48.3 x 92.4 cm).

Catalog Number	Description	Shipping Weight
80325-00	2-Basket Cart, Each basket measures 32 x 18 x 4" (81.3 x 45.7 x 10.2 cm)	43 lbs. (19.5 kg)
80450-00	4-Basket Cart, Each basket measures 16 x 18 x 4" (40.6 x 45.7 x 10.2 cm)	43 lbs. (19.5 kg)
80401-00	Extra Large Basket, 32 x 18 x 4" (81.3 x 45.7 x 10.2 cm)	5 lbs. (2.3 kg)
80402-00	Extra Small Basket, 16 x 18 x 4" (40.6 x 45.7 x 10.2 cm)	2 lbs. (0.9 kg)



PORTABLE TABLE

- Labconco's most popular cart has two shelves of phenolic board with a melamine high-pressure laminate surface which resists stains and corrosion.
- Open shelves provide easy transfer of equipment from countertop to cart.
- Dimensions: 34 7/8 x 19 x 36 3/8" high (88.6 x 48.3 x 92.4 cm).
- Shipping weight 44 lbs. (20 kg).

FLASK CARRIER CART

- Holds twenty-four 500 or 800 ml volumetric or Kjeldahl flasks either upright or neck down for easy draining.
- Top and bottom shelves of 3/8" (10 mm) fiberesin board with a melamine high-pressure laminate surface have beveled holes.
- Vinyl-coated steel neck supports hold upright flasks securely.
- Removable drip pan minimizes clean-up.
- Dimensions: 34 7/8 x 13 x 45 5/8" high (88.6 x 33 x 115.9 cm).
- Shipping weight 49 lbs. (22.2 kg).





#80210-00

NEW

REACH & ROLL CART

- Step ladder and cart in one.
- Two removable, watertight, epoxy-coated steel pans with gray vinyl mats, and one vinyl-coated wire basket.
- Sturdy, two-step ladder with non-skid surface conforms to OSHA and ANSI design guidelines.
- Spring-loaded, oversized steps with skid-resistant, rubber feet prevent cart from moving while user stands on the ladder.
- Steps withstand loads up to 1000 lbs. (453.6 kg).
- Dimensions: Cart is 19 x 39 1/2 x 43" high (48.3 x 100.3 x 109.2 cm). Top pan is 6 1/2 x 16 3/4 x 2" (16.5 x 42.5 x 5.1 cm). Bottom pan is 17 1/2 x 16 3/4 x 2" (44.5 x 42.5 x 5.1 cm). Basket is 16 x 18 x 4" (40.6 x 45.7 x 10.2 cm).
- Shipping weight 69 lbs. (31.3 kg).



#80600-00

LABORATORY BENCHES

- **Heavy Duty Auxiliary Bench Supports 540 lbs.**
- **Corrosion-Resistant Tops with Melamine High-Pressure Laminate Surfaces**
- **5" Diameter Toe-Locking Casters**

MOBILE BENCH

- Dimensions: 48 x 24 x 37 3/16" high. (121.9 x 61 x 94.5 cm). Work surface is 1 1/8" (29 mm) thick.
- Shipping weight 111 lbs. (50.3 kg).



#80750-00

VARIABLE HEIGHT BENCH

- Four adjustments from 29 5/8 to 37 3/8" (75.1 to 94.9 cm) using steel pins with nuts inserted in vertical frame members.
- Work surface dimensions: 38 x 28 x 1 1/8" thick (96.5 x 71.1 cm x 29 mm).
- Shipping weight 111 lbs. (50.3 kg).



Protecting your laboratory environment

LABCONCO



For more information, please contact us:

[ExpotechUSA](#)
[10700 Rockley Road](#)
[Houston, Texas 77099](#)
[USA](#)

[281-496-0900 \[voice\]](#)

[281-496-0400 \[fax\]](#)

E-mail: sales@expotechusa.com

Website: www.ExpotechUSA.com