

# Handheld Meters

Portable and Handheld Meters

The 315i, 330i and 340i series of handheld meters provide reliable measurements under difficult conditions both in the Lab and in the field. Lightweight and compact, the Series 3xxi meters are impact-resistant, hose-proof to IP 66 and meet the requirements for IP 67.

## Model 315i

### Simple

Simple operation means that measuring errors are avoided. All the function keys of the 315i can be operated even when wearing gloves; automatic calibration and AutoRead functions ensure stable and reproducible results. The multifunctional display for pH, oxygen, conductivity and temperature is easy to read.

## Model 330i

### Proven

The 330i, for field use with proven reliability, has a datalogger for 800 data records, multifunctional display, manual input for data, and 3000 hours of continuous operation.



## Model 340i

### All-rounder

In addition to all of the advanced features of all 3xxi series units, the model 340i can operate either with line or battery power. The model 340i also has a bi-directional interface for output to a PC or printer.



- For routine measurements
- Simple operation
- Reliable, Accurate

- Routine and standard measurements
- Wide Range of Sensors Available
- Datalogger

- Bidirectional interface RS 232
- Remote Control

## The System

### The key word is MobiLab

Our professional case can be used as a mobile laboratory. It contains a work platform with measuring beakers, a plug-in stand, all the necessary calibration and cleaning solutions as well as waterproof instructions – simply Ready to go.



# NEW

## Handheld Meter Multi 350i

- Simultaneous measurement of 3 parameters + temperature
- Back-lit graphic display
- Menu-driven

Flexible and robust. Compact precision without compromises. For simultaneous measurement of up to 4 parameters. Graphics display with user-friendly menu navigation as well as the highest level of accuracy and resolution make this instrument top class.

**Exclusively for the Multi 350i –  
the new multiparameter Probes from WTW:**

**ConOx:** combined conductivity and oxygen sensor for laboratory and field applications

**MPP 350:** outdoor multiparameter probe for pH, oxygen and conductivity

*More on page 58.*



For more information, please contact us:

[ExpotechUSA](#)  
[10700 Rockley Road](#)  
[Houston, Texas 77099](#)  
[USA](#)

[281-496-0900 \[voice\]](#)

[281-496-0400 \[fax\]](#)

E-mail: [sales@expotechusa.com](mailto:sales@expotechusa.com)

Website: [www.ExpotechUSA.com](http://www.ExpotechUSA.com)