

BOD self-check measurement – for a larger number of samples with simple sample management

OxiTop® Control 6, OxiTop® Control 12



OxiTop® Control 12

- Controller-driven
- Simultaneous measurement of up to 100 samples
- Statistical evaluation
- Automatic sample ID
- Worldwide approved method

Complete package for 6 or 12 simultaneous measurements

OxiTop® Control is the logical further development of the successful OxiTop® system using software-controlled functions and infrared interface to communicate with a handy controller, the OC 100. This system enables the simultaneous and grouped start, management, storage and tracking of 100 measuring heads via the controller and tracked on a large display with graphic evaluation. Data can be transferred to the PC for evaluation and documentation via the **AK-540/B** cable (order no. 902 842) and the communication program **Achat OC** (order no. 208 990).

Controller OC 110 in combination with the OxiTop® Control S6 / S12 is ideal for users with other applications apart from BOD (see page 76).

Check sampling progress!

The data can be called up at any time, even during sampling, thus enabling checking of the samples for errors. The display of the progress curve allows immediate detection of irregularities and interferences, such as a BOD value set too high for the volume used or undesired nitrification. Corrections can thus be made at an early stage.



For more information, please contact us:

ExpotechUSA

10700 Rockley Road
Houston, Texas 77099
USA

281-496-0900 [voice]

281-496-0400 [fax]

E-mail: sales@expotechusa.com

Website: www.ExpotechUSA.com