

Communication Program ACHAT OC for BOD Measurement with OxiTop® Control

Stored measuring data can be transferred simple and comfortable from the OxiTop® Control system to a PC with the communication program ACHAT OC. These transferred data can be processed by PC or the Controller OC 100/OC 110.


Min. software requirements:

Windows 3.1 update options on request.

Necessary accessories:

Connection cable AK 540/B.

PC Software Multi/ACHAT II: Datenkommunikation mit WTW-Photometern



- dBASE format
- Sorting according to measurement parameter
- Windows program with context-sensitive menus and help programs

WTW offers the **Multi/ACHAT II** PC software program so that stored measuring data can be transferred simply and comfortably from the instrument to a PC for documentation and statistical evaluation.

Ordering Information

MultiLab® pilot		Order No.
KOM pilot	Communications package, consisting of 1 x MultiLab® pilot and 1 AK 340/B	902 915
AK M-PC	Adapter cable for connecting a PC to the multifunction box, incl. MultiLab® pilot	902 901
AK T-PC	Interface cable for connecting the inoLab® 740/Level 3 terminal to the PC	902 902
ADA USB/Ser	Adapter USB serial interface RS 232 (9-pin socket)	902 880
Multi/ACHAT II		Order No.
Multi/ACHAT II	Software for PC under Windows, German and English	902 750
KOM Labor	Communications package, consisting of: 1 x Multi/ACHAT II and 1 AK Labor	902 754
AK Labor	Interface cable for connecting instrument to a PC via the RS 232 interface including 9/25-pin adapter and operating instructions	902 758
ACHAT OC		Order No.
ACHAT OC	PC communication program for Controller OxiTop® OC 100 or OC 110 for further processing of measuring data	208 990
AK 540/B	Interface cable for connecting instrument to a PC via the RS 232 interface including 9/25-pin adapter and operating instructions	902 842
Additional accessories see brochure "Product Details"		

For more information, please contact us:

[ExpotechUSA](#)
[10700 Rockley Road](#)
[Houston, Texas 77099](#)
[USA](#)

[281-496-0900 \[voice\]](#)

[281-496-0400 \[fax\]](#)

E-mail: sales@expotechusa.com

Website: www.ExpotechUSA.com