

DO Meters

Oxygen Sensors and Accessories

Oxygen Meters



Cellox® 325

Membrane covered galvanic oxygen sensor

- Temperature compensation
- Up to 6 months usage of membrane
- High signal resolution
- Rapid response
- Waterproof sensor (IP 68 – 2 bar)
- Incl. calibration vessel OxiCal®-SL
- Membrane leakage monitoring

Galvanic sensors: immediately ready for measurement

Equipped with protective hood for rough conditions

Simple air calibration using calibration vessel



DurOx® 325

only for handheld and field meters Oxi 3xxi/Oxi 197i and Multi 350i

Membrane covered galvanic oxygen sensor

- Temperature compensation
- Long working life – up to 6 months with one electrolyte filling
- Low approach flow
- Waterproof sensor (IP 68 – 2 bar)
- Incl. calibration vessel OxiCal®-D
- With protective hood SK-D as standard

StirrOx® G

for inoLab® Oxi 730, inoLab® Oxi 740 and ProfilLine Oxi 197i

Self-stirring oxygen sensor – simultaneous stirring and measurement

- Single-handed operation for rapid series measurements
- Constant approach flow for high reproducibility
- Extremely low self-consumption of oxygen – only 0.008 µg h⁻¹ (mg/l)⁻¹
- Incl. calibration vessel OxiCal®-ST
- Temperature compensation
- Membrane leakage monitoring

Calibration and storage vessels

Various calibration and storage vessels are available for oxygen sensors – see brochure “Product Details”.

Other sensors and accessories see brochure “Product Details”.

Ordering Information

| Oxygen sensors | (The sensor prices include the corresponding accessory case with spare parts and maintenance supplies) | Order No. |
|----------------|---|-----------|
| StirrOx® G | Self-stirring oxygen sensor for oxygen determination in Karlsruhe bottles and Winkler bottles, with OxiCal®-ST calibration and storage vessel | 201 425 |
| Cellox® 325 | Galvanic oxygen sensor with OxiCal®-SL calibration and storage vessel, waterproof plug, , cable length 1.5 m | 201 533 |
| DurOx® 325-3 | Galvanic oxygen sensor with OxiCal®-D calibration vessel, water proof plug, cable length 3 m | 201 570 |

For more information, please contact us:

[ExpotechUSA](#)
[10700 Rockley Road](#)
[Houston, Texas 77099](#)
[USA](#)

[281-496-0900 \[voice\]](#)

[281-496-0400 \[fax\]](#)

E-mail: sales@expotechusa.com

Website: www.ExpotechUSA.com